OKISTAR 6
OKISTAR hyal 6
OKISTAR 7
OKISTAR hyal 7

INSTRUCTION

**OKISTAR 6, OKISTAR hyal 6,
OKISTAR 7, OKISTAR hyal 7**
*Sterile solution for inhalation***Description**

OKISTAR 6, OKISTAR hyal 6, OKISTAR 7, OKISTAR hyal 7 sterile solutions for inhalation by nebulizer, facilitating breathing by fluidifying the secretion and improving its discharge from the mucous membrane of the upper and lower respiratory tract in cases of mucus accumulation (acute bronchiolitis, cystic fibrosis and other conditions requiring sputum dilution and evacuation from upper and lower respiratory tract) due to the osmotic effect.

Qualitative composition**OKISTAR 6**

Sodium chloride – 60 mg/ml

Water for injection

OKISTAR hyal 6

Sodium chloride – 60.0 mg/ml

Sodium Hyaluronate – 1.0 mg/ml

Water for injection

OKISTAR 7

Sodium chloride – 70 mg/ml

Water for injection

OKISTAR hyal 7

Sodium chloride – 70.0 mg/ml

Sodium Hyaluronate – 1.0 mg/ml

Water for injection

Classification of the product

Medical device – sterilized and pyrogen-free- Class IIa.

Package

4 ml in pre-filled disposable polymer containers.

Each package contains:

10 plastic containers with sterile **OKISTAR 6** or **OKISTAR hyal 6** or **OKISTAR 7** or **OKISTAR hyal 7** solution.

Or

60 plastic containers with sterile **OKISTAR 6** or **OKISTAR hyal 6** or **OKISTAR 7** or **OKISTAR hyal 7** solution.

Only the contents of the package are sterile. The medical device has been sterilized by steam sterilization.

Indications

OKISTAR 6, OKISTAR hyal 6, OKISTAR 7, OKISTAR hyal 7 are indicated for patients with bronchitis, bronchiolitis, cystic fibrosis, chronic obstructive respiratory diseases, bronchial asthma in adults and children and other conditions requiring sputum dilution and evacuation from upper and lower respiratory tract, to reduce mucosal edema, to fluidify secretion and to facilitate breathing .

Route of administration

The solution can be inhaled by means of a nebulizer using a facemask, an actuator or a nasal cannula, and following the nebulizer manufacturer's instructions for use.

Use one, 4 ml vial twice a day or according to medical prescription.

Warnings and special precautions for use

- The solution is only intended for use via inhalation. Do not inject.
- The first administration of the solution should be under the supervision of a doctor or qualified medical staff.
- In case of bronchial spasm or cough, discontinue the administration and consult a doctor.
- Do not use if the solution is not transparent and colorless.
- Do not use if the package is opened or damaged.
- The product must be used immediately after opening.
- The product is intended for single use. After use, the unused product is no longer sterile.
- Do not reuse. Once the product has been used for the first time, any residuals of the product are not suitable for the repeated use as the product is no longer sterile.
- Do not re-sterilize. Repeated sterilization may cause cross-contamination.

- Do not freeze.
- The product must not be administered orally.
- The use by children should be under the supervision of adults.
- This product includes small parts, which can pose a choking hazard. Always keep this product out of the reach of children under the age of three.

Special warnings

- Do not use after the expiry date.
- Do not use the solutions in combination with drugs or agents as this may alter their action which can impair the effectiveness of the treatment.
- If a tendency to develop dyspnea or hypersensitivity is present, OKISTAR 6, OKISTAR hyal 6, OKISTAR 7, OKISTAR hyal 7 must not be used without medical supervision.

Side effects

OKISTAR 6, OKISTAR hyal 6, OKISTAR 7, OKISTAR hyal 7 are contraindicated in case of individual hypersensitivity to any component of the product.

In persons with individual intolerance to the solution components allergic reactions are possible.

Hypertonic saline solution can cause reversible constriction of the bronchia when inhaled by very sensitive patients with bronchial asthma or a hyperreactive bronchial system. In particularly susceptible subjects, pre-medication with bronchodilators that can help prevent bronchospastic reactions is advisable at the start of the treatment. Pre-medication with bronchodilators must be carried out under the control of a doctor.

In case of the onset of bronchospasms or persistent cough, interrupt the treatment and inform the treating doctor.

How to use

For inhalation by nebulizer:

- Check the expiry date indicated on the package. Do not use after the expiry date.
- Take one of the plastic containers from the strip of 10 plastic containers.
- Open the plastic container turning the stopper.
- Use the solution in conjunction with the nebulizer following the device manufacturer's instructions for use.

Shelf life

2 years in intact packaging.

How to store

Store at a temperature from 5° to 30° C (inclusive), away from direct light and heat sources, in the adequately closed packaging.

The expiry date applies to the product when correctly stored in intact packaging.

Disposal

The product must be disposed in accordance with applicable laws on municipal waste.

Symbols used on the packaging			
	Do not reuse – for single use		Manufacturer
	Lot number		Sterilized using steam or dry heat
	Expiry date		Do not re-sterilize
	Do not use if the package is opened or damaged		Protect from light and heat
	Consult instructions for use		Temperature limits for storage
	Non - pyrogenic		Date of manufacture

Medical device according to Directive 93/42/EEC

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